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Ankle Support-Strap TA 01

Creating support in the ankle, reducing pain, reducing inflammation, creating warmth in the ankle


The ankle joint is a relatively complex joint that is prone to injuries such as sprains and dislocations, which are the cause of other injuries such as stretching and tearing of tendons and ligaments. An ankle brace is used to prevent and treat tendon strain, which by creating a suitable pressure support of several points around the joint, they create stability in the joint, which helps to treat and reduce joint pain. joint, helps to create relative stability in the joint and also reduce pain and inflammation.

Features and Benefits:

- .Pain relief by applying heat and gentle pressure
- .Pleasant pressure feeling when using

How to Use:

Pass your foot through the larger hole of the ankle strap so that the heel fits into the smaller hole. Then, cross the strap from the sole of your foot to the top of your foot and then take it in the same direction to the back of the leg and fasten it.

Ankle Support-Splint TA 11

It has two splints to support the ankle, reduce pain, and reduce inflammation

The ankle joint is a relatively complex joint that is prone to injuries such as sprains and dislocations, which are the cause of other injuries such as stretching and tearing of tendons and ligaments. Ankle splints are used to treat tendon strain. Splints around the ankle limit the range of motion of the joint. This restriction of movement creates stability in the joint. Also, splints absorb some of the forces that enter the joint. These two things heal the injury and reduce the pain.

Features and Benefits:

- .Creating two splints to support the joint


.Two splints in the ankle cause restrictions in the joint which creates stability in the wrist

.Despite the presence of ankle braces, the range of motion

.The ankle is not restricted

How to Use:

Place your foot in such a way that the heel is placed in the special hole and the splints are in line with the ankle. Then close the adhesives.

Thumb and Wrist Stabilizer - Splint Brace TW 22

Splint single splint, immobilization of the thumb, wrist support

Due to the high volume of daily activities, the ligaments and tendons of the hand get tired and eventually stretch and loosen. The thumb is also prone to injury due to its high degree of freedom. Fracture is also a common injury in the hand. One of the ways to treat these injuries is to immobilize the hand and its joints. Splint on the thumb, Thumb and Wrist Stabilizer with Splint Brace, by restricting the movement of the thumb, helps to treat stretching and loosening of the tendon and ligaments of the thumb.

Features and Benefits:

- .Limiting the bending and straightening of the thumb
- .Proper thumb support
- .Light weight material with suitable structure for anatomical fit on the hand
- .Adjustable straps for desired fit and stability
- .Limiting the bending and straightening of the thumb
- .Laminated laminated fabric with breathable fabric with 3D holes that prevents sweating.

How to Use:

Insert your thumb into the splinted area and wrap the wrist strap around your wrist.

Wrist and Palm Stabilizer-Splints TA 11
With two splints, no movement in the wrist, treatment Carpal tunnel and stretching and loosening of wrist tendons

Due to the high volume of daily activities, the ligaments and tendons of the hand


get tired and eventually stretch and relax. Fracture is also a common injury in the hand. One of the most important injuries of the wrist and forearm is carpal tunnel syndrome. This syndrome causes numbness in the tips of the fingers and their tingling. The wrist and forearm splint has an aluminum splint on the palm and a hard polymer splint on the hand. Splints on the hand and palm, by limiting the bending and straightening of the wrist joint, help to repair the stretching and laxity of the tendons and ligaments of the wrist, as well as the fusion of fractures of the wrist and palm bones.

Features and Benefits:

- .Contains an aluminum splint and a hard polymeric Splint
- .Limiting the bending and straightening of the wrist and thumb
- .Adjustable straps for desired fit and stability

How to Use:

Pass your thumb through the hole corresponding to the thumb and close the straps tightly and in such a way that the blood vessels are not under pressure.

Maintenance:

- Avoid washing the brace with hot water.
- Avoid washing the brace with washing machine.
- Wash with soap.
- Avoid pressing and ironing the brace.

Warning:

- Consult your doctor for information on the method and duration of use.
- Avoid tightening the knee brace too tightly.
- Avoid applying direct and indirect heat on it.
- Prevent joint use between people with skin disease.

OTESSY®
Products for all walks of life

User Manual

Product Ref:

TA 01 TW 22

TA 11 TW 23

Our Goal Specialized products for the care of joints and bones, while preventing damage, help in the treatment and recovery stages of the consumer's mobility and lead to maintaining health and increasing the quality of life. OTESSY, is Product for All Walks of Life.

☎ 021-4422-3675 ✉ info@otessy.com

📄 otessyofficial 🌐 www.otessy.com

Paper Num: T-R-02-119
Date of review: 1402/12/05
Review Num: 1